




HCCS Antibody

Product Code	CSB-PA010165ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53701
Immunogen	Recombinant Human Cytochrome c-type heme lyase protein (1-268AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytochrome c-type heme lyase (CCHL) (EC 4.4.1.17) (Holocytochrome c-type synthase), HCCS, CCHL
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	HCCS

Image



Western blot

All lanes: HCCS antibody at 3.54µg/ml

Lane 1: Mouse small intestine tissue

Lane 2: Mouse heart tissue

Lane 3: Mouse skeletal muscle tissue

Lane 4: HeLa whole cell lysate

Lane 5: RAW264.7 whole cell lysate

Lane 6: MCF-7 whole cell lysate

Lane 7: HepG2 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 31 kDa

Observed band size: 31 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.