




# MTX1 Antibody

<b>Product Code</b>	CSB-PA621679ESR2HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q13505
<b>Immunogen</b>	Recombinant Human Metaxin-1 protein (217-466AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Metaxin-1 (Mitochondrial outer membrane import complex protein 1), MTX1, MTX MTXN
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	MTX1

## Image



### Western blot

All lanes: MTX1 antibody at 3.78µg/ml

Lane 1: Jurkat whole cell lysate

Lane 2: A549 whole cell lysate

Lane 3: HepG2 whole cell lysate

Lane 4: MCF-7 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 52, 48, 36 kDa

Observed band size: 52, 36 kDa