




# PSTPIP1 Antibody

<b>Product Code</b>	CSB-PA018944ESR1HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O43586
<b>Immunogen</b>	Recombinant Human Proline-serine-threonine phosphatase-interacting protein 1 protein (157-416AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Proline-serine-threonine phosphatase-interacting protein 1 (PEST phosphatase-interacting protein 1) (CD2-binding protein 1) (H-PIP), PSTPIP1, CD2BP1
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	PSTPIP1

## Image



### Western blot

All lanes: PSTPIP1 antibody at 3.99 µg/ml + Raw264.7 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 48, 46 kDa

Observed band size: 48 kDa

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.