




SLC39A4 Antibody

Product Code	CSB-PA747413ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P5W5
Immunogen	Recombinant Human Zinc transporter ZIP4 protein (23-327AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc transporter ZIP4 (Solute carrier family 39 member 4) (Zrt- and Irt-like protein 4) (ZIP-4), SLC39A4, ZIP4
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	SLC39A4

Image



Western blot


All lanes: SLC39A4 antibody at 4.45 µg/ml +
Colo320 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 69, 67 kDa

Observed band size: 69, 47 kDa



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA747413ESR1HU at dilution of 1:100