




# ST3GAL5 Antibody

<b>Product Code</b>	CSB-PA890769ESR1HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9UNP4
<b>Immunogen</b>	Recombinant Human Lactosylceramide alpha-2,3-sialyltransferase protein (80-200AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Lactosylceramide alpha-2,3-sialyltransferase (EC 2.4.99.9) (CMP-NeuAc:lactosylceramide alpha-2,3-sialyltransferase) (Ganglioside GM3 synthase) (ST3Gal V) (ST3GalV) (Sialyltransferase 9), ST3GAL5, SIAT9
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Tags & Cell Markers
<b>Target Names</b>	ST3GAL5

## Image



### Western blot

All lanes: ST3GAL5 antibody at 6.2 µg/ml

Lane 1: K562 whole cell lysate

Lane 2: U251 whole cell lysate

### Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 48, 46 kDa

Observed band size: 48, 60 kDa

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.