




# STXBP2 Antibody

<b>Product Code</b>	CSB-PA618999ESR1HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q15833
<b>Immunogen</b>	Recombinant Human Syntaxin-binding protein 2 protein (1-280AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Syntaxin-binding protein 2 (Protein unc-18 homolog 2) (Unc18-2) (Protein unc-18 homolog B) (Unc-18B), STXBP2, UNC18B
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	STXBP2

## Image



### Western blot

All lanes: STXBP2 antibody at 5.88µg/ml

Lane 1: Mouse brain tissue

Lane 2: A549 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 67, 68 kDa

Observed band size: 67 kDa

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.