





TNS3 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA731561ESR1HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q68CZ2 |
| Immunogen | Recombinant Human Tensin-3 protein (1-240AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA, IHC; Recommended dilution: IHC:1:20-1:200 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | PBS with 0.02% sodium azide, 50% glycerol, pH7.3. |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Tensin-3 (Tensin-like SH2 domain-containing protein 1) (Tumor endothelial marker 6), TNS3, TEM6 TENS1 TPP |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | TNS3 |

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA731561ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA731561ESR1HU at dilution of 1:100