




CYP3A43 Antibody

Product Code	CSB-PA862059ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HB55
Immunogen	Recombinant Human Cytochrome P450 3A43 protein (244-503AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytochrome P450 3A43 (EC 1.14.14.1), CYP3A43
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	CYP3A43

Image



Western blot

All lanes: CYP3A43 antibody at 2.41µg/ml

Lane 1: Mouse liver tissue

Lane 2: Mouse skeletal muscle tissue

Lane 3: Jurkat whole cell lysate


Lane 4: 293T whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 58, 49, 29, 45 kDa

Observed band size: 49 kDa



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA862059ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.