




# INPP5K Antibody

<b>Product Code</b>	CSB-PA866277ESR1HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9BT40
<b>Immunogen</b>	Recombinant Human Inositol polyphosphate 5-phosphatase K protein (149-448AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Inositol polyphosphate 5-phosphatase K (EC 3.1.3.56) (Skeletal muscle and kidney-enriched inositol phosphatase), INPP5K, PPS SKIP
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	INPP5K

## Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA866277ESR1HU at dilution of 1:100

**Usage** For Research Use Only. Not for use in diagnostic or therapeutic procedures.