




INPP5K Antibody

Product Code	CSB-PA866277ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BT40
Immunogen	Recombinant Human Inositol polyphosphate 5-phosphatase K protein (149-448AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inositol polyphosphate 5-phosphatase K (EC 3.1.3.56) (Skeletal muscle and kidney-enriched inositol phosphatase), INPP5K, PPS SKIP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	INPP5K

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA866277ESR1HU at dilution of 1:100