




# KCNS3 Antibody

<b>Product Code</b>	CSB-PA859028ESR2HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9BQ31
<b>Immunogen</b>	Recombinant Human Potassium voltage-gated channel subfamily S member 3 protein (1-180AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Potassium voltage-gated channel subfamily S member 3 (Delayed-rectifier K(+) channel alpha subunit 3) (Voltage-gated potassium channel subunit Kv9.3), KCNS3
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Neuroscience
<b>Target Names</b>	KCNS3

## Image



### Western blot

All lanes: KCNS3 antibody at 4.97µg/ml

Lane 1: Mouse brain tissue

Lane 2: HT29 whole cell lysate

Lane 3: A549 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 56 kDa

Observed band size: 56 kDa

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.