



LIM2 Antibody

Product Code	CSB-PA012946ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55344
Immunogen	Recombinant Human Lens fiber membrane intrinsic protein (25-110AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lens fiber membrane intrinsic protein (MP18) (MP19) (MP20), LIM2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	LIM2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.