




NPHP3 Antibody

Product Code	CSB-PA768767ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z494
Immunogen	Recombinant Human Nephrocystin-3 protein (1-130AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nephrocystin-3, NPHP3, KIAA2000
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	NPHP3

Image



Western blot
 All lanes: NPHP3 antibody at 1.57 µg/ml +
 Mouse gonadal tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 151, 74, 72, 20, 27, 80, 16
 kDa
 Observed band size: 151 kDa

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.