





NPTN Antibody

Product Code	CSB-PA016028ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y639
Immunogen	Recombinant Human Neuroplastin protein (29-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neuroplastin (Stromal cell-derived receptor 1) (SDR-1), NPTN, SDFR1 SDR1
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	NPTN

Image



Western blot
 All lanes: NPTN antibody at 4.37 µg/ml + Hela whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 45, 32, 31, 38, 44 kDa
 Observed band size: 45, 50 kDa



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA016028ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.