




# NUP155 Antibody

|                            |                                                                                                  |
|----------------------------|--------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA016191ESR1HU                                                                               |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                    |
| <b>Uniprot No.</b>         | O75694                                                                                           |
| <b>Immunogen</b>           | Recombinant Human Nuclear pore complex protein Nup155 protein (1102-1391AA)                      |
| <b>Raised In</b>           | Rabbit                                                                                           |
| <b>Species Reactivity</b>  | Human                                                                                            |
| <b>Tested Applications</b> | ELISA, WB; Recommended dilution: WB:1:1000-1:5000                                                |
| <b>Form</b>                | Liquid                                                                                           |
| <b>Conjugate</b>           | Non-conjugated                                                                                   |
| <b>Storage Buffer</b>      | PBS with 0.02% sodium azide, 50% glycerol, pH7.3.                                                |
| <b>Purification Method</b> | Antigen Affinity Purified                                                                        |
| <b>Isotype</b>             | IgG                                                                                              |
| <b>Clonality</b>           | Polyclonal                                                                                       |
| <b>Alias</b>               | Nuclear pore complex protein Nup155 (155 kDa nucleoporin) (Nucleoporin Nup155), NUP155, KIAA0791 |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                                             |
| <b>Research Area</b>       | Tags & Cell Markers                                                                              |
| <b>Target Names</b>        | NUP155                                                                                           |

## Image



Western blot  
 All lanes: NUP155 antibody at 3.47µg/ml + Hela whole cell lysate  
 Secondary  
 Goat polyclonal to rabbit IgG at 1/10000 dilution  
 Predicted band size: 156, 150 kDa  
 Observed band size: 156 kDa

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.