





PLCD4 Antibody

Product Code	CSB-PA866264ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BRC7
Immunogen	Recombinant Human 1-phosphatidylinositol 4,5-bisphosphate phosphodiesterase delta-4 protein (1-290AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	1-phosphatidylinositol 4,5-bisphosphate phosphodiesterase delta-4 (hPLCD4) (EC 3.1.4.11) (Phosphoinositide phospholipase C-delta-4) (Phospholipase C-delta-4) (PLC-delta-4), PLCD4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PLCD4

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA866264ESR1HU at dilution of 1:100



Immunofluorescent analysis of A549 cells using CSB-PA866264ESR1HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.