





SLC9A6 Antibody

Product Code	CSB-PA821677ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92581
Immunogen	Recombinant Human Sodium/hydrogen exchanger 6 protein (510-669AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sodium/hydrogen exchanger 6 (Na(+)/H(+) exchanger 6) (NHE-6) (Solute carrier family 9 member 6), SLC9A6, KIAA0267 NHE6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC9A6

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA821677ESR1HU at dilution of 1:100



Immunofluorescent analysis of 293 cells using CSB-PA821677ESR1HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.