





ABCB8 Antibody

Product Code	CSB-PA868325ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NUT2
Immunogen	Recombinant Human ATP-binding cassette sub-family B member 8, mitochondrial protein (486-735AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	ATP-binding cassette sub-family B member 8, mitochondrial (Mitochondrial ATP-binding cassette 1) (M-ABC1), ABCB8, MABC1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ABCB8


Image



Western blot
 All lanes: ABCB8 antibody at 2.02 µg/ml
 + Mouse heart tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 80, 79, 76, 69, 60 kDa
 Observed band size: 80, 55 kDa



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA868325ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human bladder cancer using CSB-PA868325ESR1HU at dilution of 1:100