





# ABCB8 Antibody

<b>Product Code</b>	CSB-PA868325ESR1HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9NUT2
<b>Immunogen</b>	Recombinant Human ATP-binding cassette sub-family B member 8, mitochondrial protein (486-735AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse
<b>Tested Applications</b>	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	ATP-binding cassette sub-family B member 8, mitochondrial (Mitochondrial ATP-binding cassette 1) (M-ABC1), ABCB8, MABC1
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	ABCB8


## Image



Western blot  
 All lanes: ABCB8 antibody at 2.02 µg/ml  
 + Mouse heart tissue  
 Secondary  
 Goat polyclonal to rabbit IgG at 1/10000 dilution  
 Predicted band size: 80, 79, 76, 69, 60 kDa  
 Observed band size: 80, 55 kDa



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA868325ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human bladder cancer using CSB-PA868325ESR1HU at dilution of 1:100

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.