





# CBR4 Antibody

<b>Product Code</b>	CSB-PA850803ESR1HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q8N4T8
<b>Immunogen</b>	Recombinant Human Carbonyl reductase family member 4 protein (1-237AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Carbonyl reductase family member 4 (EC 1.-.-) (3-ketoacyl-[acyl-carrier-protein] reductase beta subunit) (KAR beta subunit) (3-oxoacyl-[acyl-carrier-protein] reductase) (EC 1.1.1.-) (Quinone reductase CBR4) (Short chain dehydrogenase/reductase family 45C member 1), CBR4, SDR45C1
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	CBR4

## Image



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA850803ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA850803ESR1HU at dilution of 1:100

