





CBR4 Antibody

Product Code	CSB-PA850803ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N4T8
Immunogen	Recombinant Human Carbonyl reductase family member 4 protein (1-237AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Carbonyl reductase family member 4 (EC 1.-.-) (3-ketoacyl-[acyl-carrier-protein] reductase beta subunit) (KAR beta subunit) (3-oxoacyl-[acyl-carrier-protein] reductase) (EC 1.1.1.-) (Quinone reductase CBR4) (Short chain dehydrogenase/reductase family 45C member 1), CBR4, SDR45C1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CBR4

Image



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA850803ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA850803ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.