





CCHCR1 Antibody

Product Code	CSB-PA004769ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TD31
Immunogen	Recombinant Human Coiled-coil alpha-helical rod protein 1 protein (483-782AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Coiled-coil alpha-helical rod protein 1 (Alpha-helical coiled-coil rod protein) (Putative gene 8 protein) (Pg8), CCHCR1, C6orf18 HCR
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	CCHCR1

Image



Western blot
All lanes: CCHCR1 antibody at 2.67 µg/ml +
293T whole cell lysate
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 89, 100, 96 kDa
Observed band size: 89 kDa



Immunohistochemistry of paraffin-embedded
human testis tissue using CSB-
PA004769ESR1HU at dilution of 1:100

