





CCHCR1 Antibody

Product Code	CSB-PA004769ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TD31
Immunogen	Recombinant Human Coiled-coil alpha-helical rod protein 1 protein (483-782AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Coiled-coil alpha-helical rod protein 1 (Alpha-helical coiled-coil rod protein) (Putative gene 8 protein) (Pg8), CCHCR1, C6orf18 HCR
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	CCHCR1


Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA004769ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA004769ESR2HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA004769ESR2HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.