





CETN1 Antibody

Product Code	CSB-PA613263DSR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12798
Immunogen	Recombinant Human Centrin-1 protein (1-172AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Centrin-1 (Caltractin isoform 2), CETN1, CEN1 CETN
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CETN1

Image



Western blot
 All lanes: CETN1 antibody at 2.7 µg/ml + Mouse gonadal tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 20 kDa
 Observed band size: 20 kDa



Immunohistochemistry of paraffin-embedded human heart tissue using CSB-PA613263DSR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.