




# DCTN6 Antibody

<b>Product Code</b>	CSB-PA006569ESR1HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O00399
<b>Immunogen</b>	Recombinant Human Dynactin subunit 6 protein (1-190AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB, IHC, IP; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200, IP:1:200-1:2000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Dynactin subunit 6 (Dynactin subunit p27) (Protein WS-3), DCTN6, WS3
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cell Biology
<b>Target Names</b>	DCTN6

## Image



### Western blot

All lanes: DCTN6 antibody at 3.5 µg/ml

Lane 1: A549 whole cell lysate


Lane 2: U251 whole cell lysate

### Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 21 kDa


Observed band size: 21, 40 kDa




### Immunohistochemistry of paraffin-embedded

human gastric cancer using CSB-

PA006569ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA006569ESR1HU at dilution of 1:100



Immunoprecipitating DCTN6 in A549 whole cell lysate

Lane 1: Rabbit control IgG instead of (1µg) instead of CSB-PA006569ESR1HU in A549 whole cell lysate. For western blotting, a HRP-conjugated light chain specific antibody was used as the Secondary antibody (1/50000)  
 Lane 2: CSB-PA006569ESR1HU (4µg) + A549 whole cell lysate (500µg)  
 Lane 3: A549 whole cell lysate (20µg)

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.