




DSN1 Antibody

Product Code	CSB-PA880969ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H410
Immunogen	Recombinant Human Kinetochores-associated protein DSN1 homolog protein (1-356AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kinetochores-associated protein DSN1 homolog, DSN1, C20orf172 MIS13
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DSN1

Image



Western blot

All lanes: DSN1 antibody at 3.09 µg/ml

Lane 1: Hela whole cell lysate

Lane 2: MCF-7 whole cell lysate

Lane 3: Mouse lung tissue


Lane 4: Mouse gonadal tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 41, 33, 39, 28 kDa

Observed band size: 41, 34, 23 kDa



Immunohistochemistry of paraffin-embedded

human brain tissue using CSB-PA880969ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.