





IL36G Antibody

Product Code	CSB-PA873706ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZH8
Immunogen	Recombinant Human Interleukin-36 gamma protein (1-169AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-36 gamma (IL-1-related protein 2) (IL-1RP2) (Interleukin-1 epsilon) (IL-1 epsilon) (Interleukin-1 family member 9) (IL-1F9) (Interleukin-1 homolog 1) (IL-1H1), IL36G, IL1E IL1F9 IL1H1 IL1RP2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL36G

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA873706ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human bladder cancer using CSB-PA873706ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.