





TTLL7 Antibody

Product Code	CSB-PA765082ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ZT98
Immunogen	Recombinant Human Tubulin polyglutamylase TTLL7 protein (1-260AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tubulin polyglutamylase TTLL7 (EC 6.-.-.-) (Testis development protein NYD-SP30) (Tubulin--tyrosine ligase-like protein 7), TTLL7
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TTLL7


Image



Western blot
 All lanes: TTLL7 antibody at 2.04µg/ml + HeLa whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 103, 78, 93 kDa
 Observed band size: 103 kDa



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA765082ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human pancreatic tissue using CSB-PA765082ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.