





RPP25L Antibody

Product Code	CSB-PA822707ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N5L8
Immunogen	Recombinant Human Ribonuclease P protein subunit p25-like protein (1-163AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribonuclease P protein subunit p25-like protein (RNase P protein subunit-like p25) (Rpp25-like protein), RPP25L, C9orf23
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RPP25L

Image



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA822707ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA822707ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.