




DHTKD1 Antibody

Product Code	CSB-PA842685ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96HY7
Immunogen	Recombinant Human Probable 2-oxoglutarate dehydrogenase E1 component DHKTD1, mitochondrial protein (1-280AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable 2-oxoglutarate dehydrogenase E1 component DHKTD1, mitochondrial (EC 1.2.4.2) (Dehydrogenase E1 and transketolase domain-containing protein 1), DHTKD1, KIAA1630
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DHTKD1

Image



Western blot

All lanes: DHTKD1 antibody at 1.81µg/ml

Lane 1: Mouse liver tissue

Lane 2: Mouse kidney tissue

Lane 3: A431 whole cell lysate

Lane 4: MCF-7 whole cell lysate


Lane 5: HepG2 whole cell lysate

Secondary


Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 103 kDa

Observed band size: 103 kDa



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA842685ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA842685ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.