




GRPEL1 Antibody

Product Code	CSB-PA875707ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HAV7
Immunogen	Recombinant Human GrpE protein homolog 1, mitochondrial protein (1-217AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	GrpE protein homolog 1, mitochondrial (HMGE) (Mt-GrpE#1), GRPEL1, GREPEL1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GRPEL1

Image



Western blot

All lanes: GRPEL1 antibody at 2.65µg/ml

Lane 1: A431 whole cell lysate


Lane 2: 293T whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 24 kDa


Observed band size: 24, 26 kDa



Immunohistochemistry of paraffin-embedded

human colon cancer using CSB-

PA875707ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA875707ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.