





MIPEP Antibody

Product Code	CSB-PA013835ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99797
Immunogen	Recombinant Human Mitochondrial intermediate peptidase protein (504-713AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitochondrial intermediate peptidase (MIP) (EC 3.4.24.59), MIPEP, MIP
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	MIPEP

Image



Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA013835ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human pancreatic cancer using CSB-PA013835ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.