





MYO5B Antibody

Product Code	CSB-PA892471ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULV0
Immunogen	Recombinant Human Unconventional myosin-Vb protein (1440-1690AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Unconventional myosin-Vb, MYO5B, KIAA1119
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MYO5B


Image



Western blot
 All lanes: MYO5B antibody at 1.7μg/ml
 Lane 1: Jurkat whole cell lysate
 Lane 2: HepG2 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 214, 112, 48 kDa
 Observed band size: 214, 112, 48 kDa



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA892471ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA892471ESR1HU at dilution of 1:100