





POP4 Antibody

Product Code	CSB-PA018373ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95707
Immunogen	Recombinant Human Ribonuclease P protein subunit p29 protein (1-220AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribonuclease P protein subunit p29 (hPOP4) (EC 3.1.26.5), POP4, RPP29
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	POP4

Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA018373ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA018373ESR1HU at dilution of 1:100