





C8B Antibody

Product Code	CSB-PA004195ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07358
Immunogen	Recombinant Human Complement component C8 beta chain protein (55-350AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Complement component C8 beta chain (Complement component 8 subunit beta), C8B
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	C8B

Image



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA004195ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human pancreatic cancer using CSB-PA004195ESR1HU at dilution of 1:100