





# FBXL5 Antibody

<b>Product Code</b>	CSB-PA892451ESR1HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9UKA1
<b>Immunogen</b>	Recombinant Human F-box/LRR-repeat protein 5 protein (1-310AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	F-box/LRR-repeat protein 5 (F-box and leucine-rich repeat protein 5) (F-box protein FBL4/FBL5) (p45SKP2-like protein), FBXL5, FBL4 FBL5 FLR1
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cell Biology
<b>Target Names</b>	FBXL5

## Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA892451ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human spleen tissue using CSB-PA892451ESR1HU at dilution of 1:100

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.