




FBXO22 Antibody

Product Code	CSB-PA854115ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NEZ5
Immunogen	Recombinant Human F-box only protein 22 protein (104-403AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	F-box only protein 22 (F-box protein FBX22p44), FBXO22, FBX22
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FBXO22

Image



Western blot

All lanes: FBXO22 antibody at 2.42 µg/ml

Lane 1: Mouse heart tissue


Lane 2: Mouse liver tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 45, 5, 31 kDa

Observed band size: 45 kDa



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-

PA854115ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.