





FBXO6 Antibody

Product Code	CSB-PA873634ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NRD1
Immunogen	Recombinant Human F-box only protein 6 protein (104-293AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	F-box only protein 6 (F-box protein that recognizes sugar chains 2) (F-box/G-domain protein 2), FBXO6, FBG2 FBS2 FBX6
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FBXO6

Image



Immunohistochemistry of paraffin-embedded human glioma using CSB-PA873634ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA873634ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.