





RNF13 Antibody

Product Code	CSB-PA019831DSR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43567
Immunogen	Recombinant Human E3 ubiquitin-protein ligase RNF13 protein (204-381AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase RNF13 (EC 2.3.2.27) (RING finger protein 13) (RING-type E3 ubiquitin transferase RNF13), RNF13, RZF
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RNF13


Image



Western blot
 All lanes: RNF13 antibody at 1.22µg/ml + MCF-7 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 43, 30 kDa
 Observed band size: 43 kDa



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA019831DSR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human glioma using CSB-PA019831DSR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.