





RNF135 Antibody

Product Code	CSB-PA811599ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IUD6
Immunogen	Recombinant Human E3 ubiquitin-protein ligase RNF135 protein (339-432AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase RNF135 (EC 2.3.2.27) (RIG-I E3 ubiquitin ligase) (REUL) (RING finger protein 135) (RING-type E3 ubiquitin transferase RNF135) (Riplet), RNF135
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RNF135

Image



Western blot
 All lanes: RNF135 antibody at 5.88µg/ml
 Lane 1: 293T whole cell lysate
 Lane 2: Jurkat whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 48, 23, 31 kDa
 Observed band size: 48 kDa



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA811599ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.