





RNF139 Antibody

Product Code	CSB-PA837849DSR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WU17
Immunogen	Recombinant Human E3 ubiquitin-protein ligase RNF139 protein (510-664AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase RNF139 (EC 2.3.2.27) (RING finger protein 139) (RING-type E3 ubiquitin transferase RNF139) (Translocation in renal carcinoma on chromosome 8 protein), RNF139, TRC8
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RNF139

Image



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA837849DSR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA837849DSR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.