





TRIM10 Antibody

Product Code	CSB-PA871394ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UDY6
Immunogen	Recombinant Human Tripartite motif-containing protein 10 protein (222-481AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tripartite motif-containing protein 10 (B30-RING finger protein) (RING finger protein 9), TRIM10, RFB30 RNF9
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TRIM10

Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA871394ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA871394ESR2HU at dilution of 1:100