




TRIM23 Antibody

Product Code	CSB-PA024459ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P36406
Immunogen	Recombinant Human E3 ubiquitin-protein ligase TRIM23 protein (1-280AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase TRIM23 (EC 2.3.2.27) (ADP-ribosylation factor domain-containing protein 1) (GTP-binding protein ARD-1) (RING finger protein 46) (RING-type E3 ubiquitin transferase TRIM23) (Tripartite motif-containing protein 23), TRIM23, ARD1 ARFD1 RNF46
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TRIM23

Image



Western blot

All lanes: TRIM23 antibody at 4.55 µg/ml +


Raw264.7 whole cell lysate

Secondary


Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 65, 64, 62 kDa

Observed band size: 65, 35 kDa



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA024459ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human thyroid tissue using CSB-PA024459ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.