





UQCR10 Antibody

Product Code	CSB-PA890670ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UDW1
Immunogen	Recombinant Human Cytochrome b-c1 complex subunit 9 protein (1-63AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytochrome b-c1 complex subunit 9 (Complex III subunit 9) (Complex III subunit X) (Cytochrome c1 non-heme 7 kDa protein) (Ubiquinol-cytochrome c reductase complex 7.2 kDa protein), UQCR10, UCRC
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	UQCR10

Image



Immunohistochemistry of paraffin-embedded human heart tissue using CSB-PA890670ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA890670ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.