




UPK3A Antibody

Product Code	CSB-PA025657ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75631
Immunogen	Recombinant Human Uroplakin-3a protein (19-207AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Uroplakin-3a (UP3a) (Uroplakin III) (UPIII), UPK3A, UPK3
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	UPK3A

Image



Western blot

All lanes: UPK3A antibody IgG at 1.76µg/ml +
 Mouse lung tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 31, 18 kDa
 Observed band size: 31 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.