





IQGAP3 Antibody

Product Code	CSB-PA771447ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86VI3
Immunogen	Recombinant Human Ras GTPase-activating-like protein IQGAP3 protein (1437-1631AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras GTPase-activating-like protein IQGAP3, IQGAP3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	IQGAP3

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA771447ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human ovarian cancer using CSB-PA771447ESR1HU at dilution of 1:100