





KLHL13 Antibody

Product Code	CSB-PA865211ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P2N7
Immunogen	Recombinant Human Kelch-like protein 13 protein (1-260AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kelch-like protein 13 (BTB and kelch domain-containing protein 2), KLHL13, BKLHD2 KIAA1309
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	KLHL13


Image



Western blot
All lanes: KLHL13 antibody at 5.06µg/ml + Jurkat whole cell lysate
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 74, 70, 73, 75 kDa
Observed band size: 74 kDa



Immunohistochemistry of paraffin-embedded human ovarian cancer using CSB-PA865211ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA865211ESR1HU at dilution of 1:100