





MCF2L Antibody

Product Code	CSB-PA521145ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15068
Immunogen	Recombinant Human Guanine nucleotide exchange factor DBS protein (850-980AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Guanine nucleotide exchange factor DBS (DBL's big sister) (MCF2-transforming sequence-like protein), MCF2L, KIAA0362 OST
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MCF2L

Image



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA521145ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human spleen tissue using CSB-PA521145ESR2HU at dilution of 1:100