





MEMO1 Antibody

Product Code	CSB-PA896500ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y316
Immunogen	Recombinant Human Protein MEMO1 protein (1-297AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein MEMO1 (C21orf19-like protein) (Hepatitis C virus NS5A-transactivated protein 7) (HCV NS5A-transactivated protein 7) (Mediator of ErbB2-driven cell motility 1) (Mediator of cell motility 1) (Memo-1), MEMO1, C2orf4 MEMO NS5ATP7
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MEMO1

Image



Immunohistochemistry of paraffin-embedded human thyroid tissue using CSB-PA896500ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA896500ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.