





RAPGEF4 Antibody

Product Code	CSB-PA820189ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WZA2
Immunogen	Recombinant Human Rap guanine nucleotide exchange factor 4 protein (1-210AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rap guanine nucleotide exchange factor 4 (Exchange factor directly activated by cAMP 2) (Exchange protein directly activated by cAMP 2) (EPAC 2) (cAMP-regulated guanine nucleotide exchange factor II) (cAMP-GEFII), RAPGEF4, CGEF2 EPAC2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RAPGEF4

Image



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA820189ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA820189ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.